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Contents

Page

1. Product description	3
1.1 Intended use of the product.	10
1.2 User groups.	10
1.3 Warranty terms	10
1.4 Declaration of conformity	10
2. Information on the safety and protection of the product	10
3. Use	11
3.1 Cautions	11
3.2 Positioning the patient on the stokes basket	11
3.3 Lifting the stokes basket	11
3.4 Lifting the stokes basket with a transport gear	12
3.5 Transport of the stokes basket on sloping surfaces or slopes, vertical transport	12
4. Cleaning and care	13
5. Maintenance	13
6. Storage	14
7. Service life	14
8. Optional accessories	14
9. Identification and warranty certificate	15

Please read the instructions for use and safety instructions carefully before you start using the ultraMINING. Keep the instructions for use in a safe place to avoid operating errors.

Should you require further information or encounter difficulties in its use or maintenance, please contact your local dealer or ultraMEDIC GmbH directly. Always quote the LOT- or serial number to identify the ultraMINING.

1.1 Intended use of the product

The ultraMINING stokes basket is a medical device and is used to safely transport of injured person.

1.2 User groups

This product may only be used by qualified, expert personnel. A prerequisite is theoretical and practical experience in handling stokes baskets.

1.3 Warranty terms

The manufacturer grants a warranty period of 24 months on this product. The warranty only covers the costs of the spare parts. Any incurred shipping or wage costs are not subject to warranty.

Please check the shipment for damage and completeness immediately upon opening. If damage is found, this must be reported within two working days, packaging damage must be reported immediately. After this period, no guarantee or warranty claims for transport damage can be asserted. No responsibility is taken for improper assembly.

1.4 Declaration of conformity

The ultraMINING was manufactured in accordance with the current EU regulation 2017/745. A declaration of conformity for these products can be acquired from the local dealership or the manufacturer.

2. Information on the safety and protection of the product

1. Before and after every use, check that the components and safety belts are intact and complete.
2. Before using the ultraMINING for the first time, conduct a test run to familiarize yourself with the product.

-
3. In the event of damage or suspected damage to the safety belts or ultraMINING remove the product from operation immediately until any occurring irregularities have been remedied.

3. Use

3.1 Cautions

1. A specialist medical assessment must be conducted before the injured party is moved, lifted or transported.
2. To ensure the patient's safety, at least 2 users are required to operate the ultraMINING.
3. The injured party must be immobilised on the ultraMINING during every use. Insufficient or no restraint may lead to severe patient or rescuer injuries.
4. Serious incidents that have occurred in connection with the product must be reported by the user and / or patient to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

3.2 Positioning the patient on the stokes basket

1. Before admitting the patient, bring the stokes basket to the best possible position, close to the patient.
2. Depending on the situation, it may be necessary to move the patient. Then transfer the injured patient from the previous immobilisation system.
3. Secure the patient to the stokes basket using the patient restraint straps. Make sure the straps are not overtightened.

3.3 Lifting the stokes basket

1. Secure the patient to the stokes basket using the patient restraint straps. Make sure the straps are not overtightened.

-
2. 2 telescopic handles are provided at the head and foot end of the stokes basket for easy lifting. To adjust the length of the handles, pull up the locking pins and pull the handles out until the locking pin sinks back into the next drill hole (image 1).


Image 1: *Length adjustment of the telescopic handles*

3.4 Lifting the stokes basket with a transport gear

1. Hang the carabiners of the transport gear in the handle openings (shorter strands at the head end) (image 2).
2. Make sure carabiners are properly locked.
3. Make sure that the transport sling is not twisted and that the injured person is securely restrained.
4. Hang the carabiners of the rescue hoist in the eye of the transport gear and slowly lift the stoker basket with the injured person.


Image 2: *Fastening of a transport gear*

3.5 Transport of the stokes basket on sloping surfaces or slopes, vertical transport

1. Secure the patient with at least 4 straps, 2 of them under the armpits, to prevent them from slipping. If the patient is injured in the lower limbs, immobilize them with suitable

additional equipment.

2. Make sure the straps are not overtightened.
3. Before lifting or lowering the stokes basket, attach it to a static safety system.

ATTENTION! Please note the instructions for attaching the transport gear. As an optional accessory, we recommend ultraCABLE (SAN-0099) from ultraMEDIC. Use of a transport gear from other manufacturers is at your own risk.

4. Cleaning and care

Use warm water (40°C) to clean the ultraMINIG, the straps and the rescue sheet.

For disinfection, you can use commercially available, alcohol-based agents. Other solutions or cleaning chemicals may damage the material of the ultraMINING and straps.

Cleaning, care and disinfecting agents that contain chlorine, solid abrasives, acids or alkalis should not be used, as all of these substances can attack and destroy the protective oxide layer of the stokes basket.

5. Maintenance

In order to ensure the safety of the injured party and the user, we recommend that you have the ultraMINING visually inspected by an expert at least **once a year**.

For users in Germany, we recommend consideration of the following regulations and guidelines:

DGUV Grundsatz 305-002

DGUV Regel 112-198/199

MPBetreibV

The following parts must be checked and the following defects must be excluded:

Base plate: deformation, tears, breaks, loosened screws, broken parts

Closures: impaired function, deformation, torsion, tears and breaks

Straps: Cutting and abrasion points, damaged seams, tissue-damaging stains

Metal parts: Rust and oxidation

other damages

In the event of damage or suspected damage, remove the product from operation immediately.

7. Storage

Store the ultraMINING in a dry and dark place. Avoid direct sunlight and direct heat radiation from heat sources. These factors could damage the material.

8. Service life

The discard criteria depends on the product, its frequency of use and the external conditions of use. The service life of the ultraMINING depends on the individual conditions of use. Depending on the frequency and intensity of use of the transportation suspension, certain wear and tear is normal. Appropriate use is assumed here.

With regard to the straps, we recommend replacing it after 6 years at the latest.

The service life of the optional accessories can be found in the associated instructions for use.

9. Optional accessories

SAN-0260-2	X-BOON TWO	Vacuum mattress
SAN-0278	ultraVACPUMP I	Vacuum pump
SAN-0099	ultraCABLE	Transport gear made of stainless steel
SAN-0157		Rescue sheet
BAUG-143		Replacement safety belts, set

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