

⊕ ultraTOWER



ultraTOWER

Stationäre Säule zur präventiven Händedesinfektion

- extrem stabil und standfest
- aus pulverbeschichtetem Metall
- geeignet für 500 und 1000 ml Flaschen
- Flaschensicherung
- vollautomatischer kontaktloser Sensor-Spender mit Auffangschale
- in 4 Stufen verstellbare Sprühmenge und Sprühdauer
- eingefrästes Logo*

Bestell-Nr.: SAN-0050
(Aseptoman Plus ist im Lieferumfang nicht enthalten)

Aseptoman Plus

Hygienische und chirurgische Händedesinfektion (EN 1500 und EN 12971)

- bakterizid inkl. MRSA
- viruzid gem. RKI-Empfehlung und EN 14476
- Noro-Viren wirksam bereits nach 15 Sek.
- parfümfrei
- rückfettend
- biozide Händedesinfektion
(fällt nicht unter das Arzneimittelgesetz, es wird keine AM-Großhandelserlaubnis benötigt)
- hygienische Händedesinfektion gem. VAH / DGHM / EN1500 konz. 30 Sek.



Bestell-Nr.:	SAN-0051	VPE: 1 x 500 ml
	SAN-0051-20	VPE: 20 x 500 ml
	SAN-0052	VPE: 1 x 1000 ml
	SAN-0052-10	VPE: 10 x 1000 ml

* Ab einer Bestellmenge von 10 Stk. kann Ihr Logo kostenneutral in die Säule eingefräst werden.

⊕ ultraTOWER



ultraTOWER

Stationary column for preventive hand disinfection

- Extremely stable
- Made of powder-coated metal
- Suitable for 500 und 1000 ml bottles
- Bottles securing
- Fully-automatic, contactless sensor-dispenser with a drip tray
- Spraying quantity and spraying time can be dispensed in 4 stages
- Milled logo*

Order-No.: SAN-0050
(Aseptoman Plus is not included in the scope of delivery)

Aseptoman Plus

Hygienic and surgical hand disinfection (EN 1500 und EN 12971)

- Bactericidal incl. MRSA
- Virucidal in accordance with RKI-Recommendation EN 14476
- Effective against Noro-viruses in 15 seconds
- Fragrance-free
- Refatting effect
- biocide hand disinfection
(falls not under the German pharmaceutical law)
- hygienic hand disinfection in accordance with VAH / DGHM / EN1500 conc. 30 sec.



Order-No.:	SAN-0051	PU: 1 x 500 ml
	SAN-0051-20	PU: 20 x 500 ml
	SAN-0052	PU: 1 x 1000 ml
	SAN-0052-10	PU: 10 x 1000 ml

* From an order quantity of 10 pieces: your company logo can be milled on the column (cost-neutral).

